

IMPORTANT INFORMATION TO READ and RETURN

Installation Requirements for a MEDIAWEL 30

Thank you for choosing one of our products for your laboratory. To enable our engineers to perform an efficient, trouble-free installation please study, complete and email this form to us on **service@dwscientific.co.uk**. Should you have any questions, please do not hesitate to contact us, as we are here to help. When we have received the completed form, our Service Department will contact you to arrange a mutually convenient installation date.

**The following information represents the ideal requirement.
Please contact us IMMEDIATELY if your intended location does not match this specification.**



<p>Access Requirements</p> <p>The unit can fit through a door opening of 700mm (27.5") but check size of lifts, stairs etc.</p> <div style="text-align: right;"><input type="checkbox"/></div>																							
<p>Space Requirements</p> <p>The weight of the equipment is 150kg.</p> <p>The unit should be sited as close as possible (preferably to the right) of the pourer stacker</p> <p>External Dimensions</p> <table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th>Width</th> <th>Depth</th> <th>Height</th> </tr> <tr> <th><i>mm</i></th> <th><i>mm</i></th> <th><i>mm</i></th> </tr> </thead> <tbody> <tr> <td>585</td> <td>735</td> <td>1053</td> </tr> </tbody> </table> <div style="text-align: right;"><input type="checkbox"/></div>			Width	Depth	Height	<i>mm</i>	<i>mm</i>	<i>mm</i>	585	735	1053												
Width	Depth	Height																					
<i>mm</i>	<i>mm</i>	<i>mm</i>																					
585	735	1053																					
<p>Mains Requirements</p> <table style="width: 100%;"> <tr> <td style="width: 50%;">Electricity Supply</td> <td style="width: 50%;">Socket</td> </tr> <tr> <td>400 +/-10% V, three- phase neutral & earth</td> <td>Five pin 16 amp</td> </tr> </table> <div style="text-align: right;"><input type="checkbox"/></div>			Electricity Supply	Socket	400 +/-10% V, three- phase neutral & earth	Five pin 16 amp																	
Electricity Supply	Socket																						
400 +/-10% V, three- phase neutral & earth	Five pin 16 amp																						
<p>Water System Requirements</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 30%;">Input (Tap water)</th> <th style="width: 30%;"></th> <th style="width: 40%;">Output</th> </tr> </thead> <tbody> <tr> <td>Minimum flow</td> <td>5 l/min</td> <td>Copper or PVC pipe that may hold high temperatures (130°C)</td> </tr> <tr> <td>Pressure</td> <td>3-4 bar</td> <td>Pipe minimum diameter: 40mm</td> </tr> <tr> <td>Temperature</td> <td>5 - 20°C</td> <td>No pressure at output</td> </tr> <tr> <td>Maximum water hardness</td> <td>16°TH</td> <td></td> </tr> <tr> <td>Water tap location</td> <td></td> <td>Located less than 2 metres from the place of usage</td> </tr> <tr> <td>Water evacuation point</td> <td></td> <td>Located less than 2 metres from the place of usage and prepared to resist temperatures up to 120°C (copper or CPVC). Is this evacuation equipped with a U-bend / siphon?</td> </tr> </tbody> </table> <div style="text-align: right;"><input type="checkbox"/></div>			Input (Tap water)		Output	Minimum flow	5 l/min	Copper or PVC pipe that may hold high temperatures (130°C)	Pressure	3-4 bar	Pipe minimum diameter: 40mm	Temperature	5 - 20°C	No pressure at output	Maximum water hardness	16°TH		Water tap location		Located less than 2 metres from the place of usage	Water evacuation point		Located less than 2 metres from the place of usage and prepared to resist temperatures up to 120°C (copper or CPVC). Is this evacuation equipped with a U-bend / siphon?
Input (Tap water)		Output																					
Minimum flow	5 l/min	Copper or PVC pipe that may hold high temperatures (130°C)																					
Pressure	3-4 bar	Pipe minimum diameter: 40mm																					
Temperature	5 - 20°C	No pressure at output																					
Maximum water hardness	16°TH																						
Water tap location		Located less than 2 metres from the place of usage																					
Water evacuation point		Located less than 2 metres from the place of usage and prepared to resist temperatures up to 120°C (copper or CPVC). Is this evacuation equipped with a U-bend / siphon?																					
<p>Connection</p> <p>¾ " bsp male (washer connection) fitting or a ½ " (15mm) barb to receive a flexible tubing (12.7mm inner diameter). In both cases, an isolation tap is recommended.</p> <div style="text-align: right;"><input type="checkbox"/></div>																							
<p>Drain</p> <p>1 ¼ " (30mm) or 1 ½ " (38mm) waste trap or a sufficiently large bore waste pipe to accommodate one 1" (25mm) hose.</p> <div style="text-align: right;"><input type="checkbox"/></div>																							

Validation is highly recommended for each load to be utilised in order to:

- ensure that the required temperatures and times throughout the load are achieved
and
- determine cycle settings.

As Don Whitley Scientific is **UKAS accredited**, we can offer you validation and calibration of your heat sterilisation equipment.

If you would like to discuss this, please call us on 01274 595728.

Delivery and installation are free of charge (unless otherwise agreed). If our engineers are unable to install the unit and a return journey is necessary, **a charge may be made.**

It is essential that this form is completed and returned to avoid delay to your installation.

THANK YOU FOR THINKING WHITLEY

Signature

Title

Print Name

Establishment